

## Hemoglobin Variants and Derivatives on A1C Measurement

Abnormal hemoglobin (Hb) variants (i.e. HbS, HbC, and HbE) and chemically modified derivatives (i.e. carbamyl-Hb) can cause interference and dramatically affect the results of A1C measurement. Such problems should be suspected whenever A1C results are out of keeping with results of self-monitoring of blood glucose. The Cholestech GDX™ System has virtually no interference from hemoglobin variants because the method uses boronate affinity chromatography to separate the glycated hemoglobin fraction from the non-glycated fraction.

Nearly 8% of African Americans carry the HbS trait and 2.3% carry HbC. In sub-Saharan Africa, prevalence of these two is up to one-third of all patients. HbE can be as high as 30% in Southeast Asia. HbF can reach 30% in individuals with hereditary persistence. Chemically modified Hbs may be chronically present in diabetic patients. Carbamylated Hb is the most commonly encountered of these.<sup>2</sup>

Because different laboratories use different methods, practitioners should familiarize themselves with the methodology of their local laboratory. Preferably a single laboratory should be used for all measurements for patients of a given physician to lessen confusion. The table below references some of the more common methods and the Hb variants that will cause interference with the measurements.

Methods	Interference (Yes/No)				
	Carbamyl-Hb	Hb C trait	Hb E trait	Hb S trait	Elevated HbF
<b>Boronate Affinity Chromatography</b>					
<b>Cholestech GDX™<sup>1</sup></b>	<b>No</b>	<b>No</b>	<b>No</b>	<b>No</b>	<b>No</b>
Primus CLC 330™ <sup>2,3</sup>	No	No	No	No	No
Abbott IMx® <sup>2,3</sup>	No	No	--	No	No
<b>Cation-exchange chromatography</b>					
Bio-Rad Diamat® <sup>2,3</sup>	Yes	Yes	No	Yes	Yes
Bio-Rad Variant® <sup>2,3</sup>	Yes	No	No	Yes	Yes
Menarini 8140 <sup>2,4</sup>	Yes	No	Yes	Yes	No
Tosoh A1c 2.2 Plus™ <sup>2</sup>	No	No	Yes	No	No
<b>Immunoassay</b>					
Bayer DCA 2000® <sup>2,3</sup>	No	No	No	No	Yes
Roche TinaQuant II® <sup>2,3</sup>	No	No	No	No	No
Roche Unimate® <sup>2</sup>	No	No	--	Yes	No

-- No Data Available

### References

1. John WG. Glycated haemoglobin analysis. *Ann Clin Biochem* 1997; 34 ( Pt 1):17-31.
2. Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycohemoglobin. *Clin Chem* 2001; 47:153-63.
3. Goldstein DE, Little RR, Lorenz RA, Malone JI et al. Tests of glycemia in diabetes. *Diabetes Care* 1995; 18:896-909.
4. Roberts WL, De BK, Brown D, Hanbury CM et al. Effects of hemoglobin C and S traits on eight glycohemoglobin methods. *Clin Chem* 2002; 48:383-5.

Cholestech GDX is a trademark of Cholestech Corporation.

©2002 Cholestech Corporation • 3347 Investment Blvd., Hayward, CA 94545-3808 USA • Phone 800-733-0404 • Fax 510-732-7227 • www.cholesteck.com  
MKT 12421 Rev. A