


# INTRODUCING THE "NXSTEP" DRUGCHECK® TEST CUP



JUST COLLECT THE SPECIMEN  
RIGHT INTO THE CUP



PEEL THE PRIVACY LABEL



READ RESULTS IN JUST 5  
MINUTES.

## Accuracy:

Our FDA-approved product has the same degree of accuracy when compared to SAMHSA Certified Laboratory screening methods. Independent studies prove test strips do not contaminate the specimen.

## Privacy Label:

DrugCheck test cup screens for adulterants using the following parameters: pH, Specific Gravity, Nitrate and Creatinine. A pull-tab label conceals test results, so only appropriate personnel can read them, adding privacy during testing.

## Self Contained:

A positive specimen can be shipped in the same container to a laboratory for GC/MS confirmation. Each cup has an 18-month shelf-life and includes an affixed temperature strip, identification label, and security seal.

AMP	Amphetamine	1000 ng/ml
BAR	Barbiturate	300 ng/ml
BZO	Benzodiazepine	300 ng/ml
COC	Cocaine	300 ng/ml
THC	Marijuana	50 ng/ml
MTD	Methadone	300 ng/ml
MET	Methamphetamine	1000 ng/ml
OPI300	Opiates 300	2000 ng/ml
OPI2000	Opiates 2000	2000 ng/ml
PCP	Phencyclidine	25 ng/ml
OXY	Oxycodone	100 ng/ml
PPX	Propoxyphene	300 ng/ml
MDMA	Ecstasy	500 ng/ml
TCA	Tricyclic Antidepressants	1000 ng/ml

- ▶ No Assembly Required
- ▶ No Moving Parts to Use
- ▶ No Dipping Of Cards
- ▶ No Tilting To Activate
- ▶ No Split Samples Or Chambers To Rely On
- ▶ No Cross Contamination to Worry About

Health Management Systems, Corp.  
6913 K. Avenue, Suite #304  
Plano, TX 75074  
888-300-4672